

Human GAGE2A Recombinant Protein

Catalog Number: PROTQ6NT46

Overview

Product Name	Human GAGE2A Recombinant Protein
Description	Human GAGE2A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-116aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Liquid. In Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q6NT46, 1-116aa; MSWRGRS TYRPRRRYV EPEMIGPMR PEQFSDEV EPATPEEGEPAT QRQDPAAAQE GQDEGASAGQ
GPKPEAHSQE QGHPQTGCEC EDGPDGQEMD PPNPEEVKTP EEGEKQSQC

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GAGE2A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.