

Mouse Sepsecs Recombinant Protein

Catalog Number: PROTQ6P6M7

Overview

| | |
|--------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name | Mouse Sepsecs Recombinant Protein |
| Description | Mouse Sepsecs Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-504aa. Application(s): SDS-PAGE. |
| Size | Starting from 4ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 0.2M NaCl, 50% glycerol |

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q6P6M7, 1-504aa; MGSSHHHHHH SSGLVPRGSH MGSMNPESFA AGERRVSPAY VRQGCEARRA HEHLIRLLLE QGKCPEDGWD
ESTLEFLHE LAVMDSNNFL GNCVGEREG RVASALVARR HYRFIHGIGR SGDISAVQPK AAGSLLNKI TNSLVLNVIK LAGVHSVASC
FVVPMTGMS LTLCFLTRH KRPKAKYIHW PRIDQKSCFK SMVTAGFEPV VIENVLEGDE LRTDLKAVEA KIQELGPEHI LCLHSTTACF
APRVPDRLEE LAVICANYDI PHVVNNAYGL QSSKCMHLIQ QGARVGRIDA FVQSLDKNFM VPVGGAIAG FNEPFIQDIS KMYPGRASAS
PSLDVLITLL SLGCSGYRKL LKERKEMFVY LSTQLKLAE AHNERLLQTP HNPISLAMTL KTIDGHHDKA VTQLGSMFLFT RQVSGARAVP
LGNVQTVSGH TFRGFMHAD NYPCAYLNAA AAIGMKMQDV DLFIKRLDKC LNIVRKEQTR ASVVSADRN KAEDADIEEM ALKLDDVLGD
VGQGPAL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Sepsecs Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.