

Human GIMAP6 Recombinant Protein

Catalog Number: PROTQ6P9H5

Overview

Product Name	Human GIMAP6 Recombinant Protein
Description	Human GIMAP6 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-292aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 30% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q6P9H5, 1-292aa; MGSSHHHHHH SGLVPRGSH MGSMEEEYE QIPQENPPEE LSQDPVLELS GGLREKEQKT PRRLRLILMG
KTGSGKSATG NSILGRDVF SKLSTRPVTK TSQRRSREWA GKELEVIDTP NILSPQVSPE VADAICQAI VLSAPGPHAVL LVTQLGRFTD
EDQQVRRRLQ EVFGVGV LGH TILVFTRKED LAGGSLEDYV RETNNQALAW LDVTLARRHC GFNNRAQGEE QEAQLRELME KVEAIMWENE
GDYYSNKAYQ YTQQNFRLKE LQERQVSQGG GSEDVPGES WLEGLS QIQK EEEAHRCLL GKADL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GIMAP6 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.