

Human Methionine Aminopeptidase 1D Recombinant Protein

Catalog Number: PROTQ6UB28

Overview

Product Name	Human Methionine Aminopeptidase 1D Recombinant Protein
Description	Human Methionine Aminopeptidase 1D Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 20-335aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 20% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q6UB28, 20-335aa; SSPLNHI YLHKQSSSQ RRNFFRRQR DISHSIVLPA AVSSAHPVPK HIKKPDYVTT GIVPDWGDSEVKNEDQIQGLHQACQLARH VLLLAGKSLK VDMTTEEIDA LVHREIISHN AYPSPLYGG FPKSVCTSVN NVLCHGIPDS RPLQDGDIIIN IDVTVYNYG YHGDTSETFLV GNVDECGKKL VEVARRCRDE AIAACRAGAP FSVIGNTISH ITHQNGFQVC PHFVGHGIGS YFHGHPEIWH HANDSDLPME EGMFTIEPI ITEGSPEFKV LEDAWTVVSL DNQRSAQFEH TVLITSRGAQ ILTKLPHEA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Methionine Aminopeptidase 1D Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.