

# CYP2R1 (NM\_024514) Human Recombinant Protein

Catalog Number: PROTQ6VVX0

## Overview

|              |   |
|--------------|---|
| Product Name | CYP2R1 (NM_024514) Human Recombinant Protein  |
| Description  | Recombinant protein of human cytochrome P450, family 2, subfamily R, polypeptide 1 (CYP2R1) |
| Size         | 20 µg   |
| Tag          | C-Myc/DDK   |
| Form         | Frozen Solution in PBS Buffer   |
| Source       | HEK293T   |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol  |

## Concentration

>50 ug/mL as determined by microplate BCA method

## Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

## Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

## Amino Acid Sequence

MWKLWRAEEGAAALGGALFLLLALGVRLKQRRPMGFPPGPPGLPFIGNIYSLAASSELPHVYMRKQSQVYGEIFSLDLGGISTVVLNGYDVVKECLV  
HQSEIFADRPCLPLFMKMTKMGLLNSRYGRGWVDHRR LAVNSFRYFGYGQKS FESKILEETKFFND AIETYKGRPFDFKQLITNAVSNITNLIIFGERFT  
YEDTDFQHMIELFSENVELAASASVFLYNAPFWIGILPFGKHQQLFRNAAVVYDFLSRLIEKASVNRKPQLPQH FVDAYLDEMDQGKNDPSSTFSKENLI  
FSVGELI IAGTETTTNVLRWAILFMALYPNIQGQVQKEIDLIMGPNGKPSWDDKCKMPYTEAVLHEVLRFCNIVPLGIFHATSEDAVVRGYSIPKGTTVIT  
NLYSVHFDEKYWRDPEVFHPERFLDSSGYFAKKEALVPFSLGRRHCLGEHLAR MEMFLFFTALLQR FHLHFPHELVPDLKPRLGMTLQPQPYLICAERR  
?/description\_after\_attributes?

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



CYP2R1 (NM\_024514) Human Recombinant Protein