

## Human recombinant Carbohydrate sulfotransferase 3/CHST3 protein

Catalog Number: PROTQ7LGC8

### Overview

Product Name	Human recombinant Carbohydrate sulfotransferase 3/CHST3 protein
Description	Human recombinant Carbohydrate sulfotransferase 3/CHST3 protein expressed in Baculovirus with His-Tag. Sequence domain: 39-479aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q7LGC8, 39-479aa; EKENKII SRVSDKCLKQI PQALADANST DPALILAENA SLLSLELDS AFSQLQSRLR NLSLQLGVPE AMEAAGEEEE  
EQRKEEPPR PAVAGPRRHV LLMATTRTGS SFVGEFFNQQ GNIFYLFEPL WHIERTVSFE PGGANAAGSA LVYRDVLKQL FLCDLVLEH  
FITPLPEDHL TQFMFRRGSS RSLCEDPVCT PFVKKVFEKY HCKNRRCGPL NVTLAAEACR RKEHMALKAV RIRQLEFLQP LAEDPRDLR  
VIQLVRDPRA VLASRMVAFA GKYKTWKKWL DDEGQDGLRE EEVQRLRGNC ESIRLSAELG LRQPAWLRGR YMLVRYEDVA RGPLQKAREM  
YRFAGIPLTP QVEDWIQKNT QAAHDGSGIY STQKNSSEQF EKWRFSMPFK LAQVVQAACG PAMRLFGYKL ARDAAALTNR SVSLLLEERT FWVT

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant Carbohydrate sulfotransferase 3/CHST3 protein

For Research Use Only. Not for use in diagnostic procedures.