

Dymeclin (DYM) (NM_017653) Human Recombinant Protein

Catalog Number: PROTQ7RTS9

Overview

| | |
|--------------|--|
| Product Name | Dymeclin (DYM) (NM_017653) Human Recombinant Protein |
| Description | Recombinant protein of human dymeclin (DYM) |
| Size | 20 µg |
| Tag | C-Myc/DDK |
| Form | Frozen Solution in PBS Buffer |
| Source | HEK293T |
| Formulation | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol |

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MGSNSSRIGDLPKNEYLKKLSGTESISENDPFWNQLLSFSFPAPTSSSELKLEEATISVCRSLVENNPRTGNL GALIKVFLSRTKELKLSAECQNHIFIWQTH
NALFIICLLKVFICQMSEELQLHFTYEEKSPGNYSSDSEDLLELLCCLMQLITDIPLLDITYEISVEAISTMVVFLSCQLFHKEVLRQSISHKYLMRGPCLPY
TSKLVKTLLYNFIRQEKP PPGAHVFPQQSDGGGLLYGLASGVATGLWTVFTLGGVGSKAAASPELSSPLANQSLLLLLVLANLTDASDAPNPYRQAIMSF
KNTQDSSPFSSIPHAFQINFNSLYTALCEQQTSDQATLLLYTLLHQNSNIRTYMLARTDMENLVLPILYIYHVEERNSSHVYMALIILLITEDYGFNRSIH
EVILKNITWYSERVLTEISLGSLLILVVIRTIQYNMTRTRDKYLHTNCLAALANMSAQFRSLHQYAAQRIISLFSLLSKKHNVLEQATQSLRGSLSNDVPLP
DYAQDLNVIEEVRMMLEIINSCLTNSLHHNPNLVYALLYKRDLEQFRTHPSFQDIMQNI DLVISFFSSRLLQAGAELSVERVLEI IKQGVVALPKDRLKKF
PELKFKYVEEQPEEFIPYVWSLVYNSAVGLYWNPQDIQLFTMDS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Dymeclin (DYM) (NM_017653) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.