

Human MAVS Recombinant Protein

Catalog Number: PROTQ7Z434

Overview

Product Name	Human MAVS Recombinant Protein
Description	Human MAVS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-513aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q7Z434, 1-513aa; MPFAEDK TYKYICRNFS NFCNVDVVEI LPYLPCLTAR DQDRLRATCT LSGNRDTLWH LFNTLQRRPG WVEYFIAALR
GCVELDLADE VASVYQSYQP RTSDRPPDPL EPPSLPAERP GPPTPAAAHS IPYNCREKE PSYPMPVQET QAPESPGENS EQALQTLSPR
AIPRNPDGGP LESSDLAAL SPLTSSGHQE QDELGSTHT AGATSSLTPS RGPVSPVSF QPLARSTPRA SRLPGPTGSV VSTGTSFSS SPGLASAGAA
EGKQGAESDQ AEPIICSSGA EAPANLPSK VPTTLMVNT VALKVPANPA SVSTVPSKLP TSSKPPGAVP SNALNPAPS KLPINSTRAG
MVPSKVPTSM VLTKVASTV PTDGSSRNEE TPAAPT PAGA TGGSSAWLDS SSENRLGSE LSKPGVLASQ VDSPFSGCFE DLAISASTSL
GMGPCHGPEE NEYKSEGTFG IHVAENPSIQ LLEGNPGPPA DPDGGPRPQA DRKFQEREVP CHRPS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MAVS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.