

Human LYG2 Recombinant Protein

Catalog Number: PROTQ86SG7

Overview

| | |
|--------------|--|
| Product Name | Human LYG2 Recombinant Protein |
| Description | Human LYG2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 20-212aa. Application(s): SDS-PAGE. |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% glycerol, 1mM DTT |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q86SG7, 20-212aa; MGSSHHHHHH SGLVPRGSH MGSSYPFSHS MKPHLHPRLY HGCYGDIMTM KTSGATCDAN SVMNCGIRGS
EMFAEMDLRA IKPYQTLIKE VGQRHCVDPA VIAAISRES HGGSVLQDGW DHRGLKFGLM QLDKQTYHPV GAWDSKEHLS QATGILTERI
KAIQKKFPTW SVAQHLKGGL SAFKSGIEAI ATPSIDNDF VNDIIRAKF YKRQSF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human LYG2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.