

## Human PHF13 Recombinant Protein

Catalog Number: PROTQ86Y18

### Overview

Product Name	Human PHF13 Recombinant Protein
Description	Human PHF13 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-300aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2mM DTT, 1mM EDTA, 250mM Imidazole

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q86Y18, 1-300aa; MGSSHHHHHH SGLVPRGSH MGSMDSDESCA AAFHPPEEYSP SCKRRRTVED FNKCTFVLA YAGYIYPYKELPLRSPSP ANSTAGTIDS DGWDAGFSDI ASSVPLPVSD RCFSHLQPTL LQRAKPSNFL LDRKKTDLK KKKKRKRKRRDS DAPGKEGYRGLLLKLEAADP YVETPTSPTL QDIPQAPSDP CSGWSDTTPS SGSCATVSPD QVKEIKTEGK RTIVRQGKQV VFRDEDSTGN DEDIMVDSDD DSWDLVTCFC MKPFAGRPMI ECNECHTWHI LSCAKIRKSN VPEVFVCQKC RDSKFDIRRS NRSRTGSRKL FLD

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human PHF13 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.