

Human HOMER1 Recombinant Protein

Catalog Number: PROTQ86YM7

Overview

Product Name	Human HOMER1 Recombinant Protein
Description	Human HOMER1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-354aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 30% glycerol, and 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q86YM7, 1-354aa; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGS HMGE QPIFSTRAHV FQIDPNTKKN WVPTSKHAVT
VSYFYDSTRN VYRIISLDGS KAIINSTITP NMTFTKTSQK FGQWADSRAN TVYGLGFSSE HHLSKFAEFK QEFKEAARLA KEKSQEKMEL
TSTPSQESAG GDLQSP LTP E SINGTDDERT PDVTQNSEPR AEPTQNALPF SHSSAISKHW EAELATLKG N NAKLTAALLE STANVKQWKQ
QLAAYQEEAE RLHKRVTELE CVSSQANAVH THKTELNQTI QELEETLKLK EEEIERLKQE IDNARELQEQ RDSL TQKLQE VEIRNKDLEG
QLSDLEQRLE K SQNEQEA FR NNLKTLLEIL DGKIFELTEL RDNLAKLLEC S

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HOMER1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.