

Mouse ALT2/GPT2 Recombinant Protein

Catalog Number: PROTQ8BGT5

Overview

Product Name	Mouse ALT2/GPT2 Recombinant Protein
Description	Mouse ALT2/GPT2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-522aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 50ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 20% glycerol, 2mM DTT

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Mouse, Q8BGT5, 1-522aa; QRAAVLVRR GSCPRASGPW GRSHSSAAAE ASAALKVRPE RSPRDRILTL ESMNPQVKAV EYAVRGPIVL KAGEIEMELQ
RGIKKPFTEV IRANIGDAHA MGQQPITFLR QVMALCTYPN LLNSPSFPED AKKRARRILQ ACGGNSLGSY SASQGVNCIR EDVAAFITRR
DGVPADPDNI YLTTGASDGI STILKLLVSG GGKSR TGVM I PIPQYPL YSA VISELDAVQV NYLDEENCW ALNVDEL RRA LRQAKDHCDP
KVLCIINPGN PTGQVQSRKC IEDVIHFAWE EKLFL LADEV YQDNVYSPDC RFHSFKKVL Y QMGHEYSSNV ELASFHSTSK GYMGECEGYRG
GYMEVINLHP EIKGQLVKLL SVRLCPPVSG QAAMDIVVNP PEPGEESFEQ FSREKEFVLG NLAKKAKLTE DLFNQVPGIQ CNPLQGAMYA
FPRIIPAKA VEAAQSHKMA PDMFYCMKLL EETGICV VPG SGFGQREGTY HFRMTILPPV DKLKTVLHKV KDFHLKFLEQ YS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse ALT2/GPT2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.