

## Mouse Ces2e Recombinant Protein

Catalog Number: PROTQ8BK48

### Overview

Product Name	Mouse Ces2e Recombinant Protein
Description	Mouse Ces2e Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 27-559aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q8BK48, 27-559aa; QDSASPIRNT HTGQVRGSLV HVKDTDIAVH TFLGIPFAKP PVGPLRFAPP EAPEPWSGVR DGTSHPNMCL QNDNLMGSED LKMMNLILPP ISMEDCLYL NIYVPAHAHE GSNLPVMVWI HGGALTVGMA SMYDGSMMLAA TEDVVVVAIQ YRLGVLGFFS TGDQHAKGNW GYLDQVAALR WVQQNIVHFG GNPDRVTIFG ESAGGTSVSS HVVSPMSQGL FHGAIMESGV AVLPDLISS SEMVHRIVAN LSGCAAVNSE TLMCCLRGKN EAEMLAINKV FKII PGVVDG EFLPKHPQEL MASKDFHPVP SIIGINNDEY GWILPTIMDP AQKIEEITRK TLPVVKSTA LKMMPLPECG DLLMEEYMGD TEDPETLQAQ FREMKGDFMF VIPALQVAHF QRSHAPVYFY EFQHRPSFFK DFRPPYVKAD HGDEIFLVFG YQFGNIKLPY TEEEEQLSRR IMKYWANFAR HGNPNSEGLP YWPVMDHDEQ YLQLDIQPSV GRALKARRLQ FWTKTLPQKI QELKGSQERH KEL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Ces2e Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.