

Mouse Galectin-3 Recombinant Protein

Catalog Number: PROTQ8C253

Overview

Product Name	Mouse Galectin-3 Recombinant Protein
Description	Mouse Galectin-3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-264aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% glycerol, 1mM DTT, 2mM EDTA

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q8C253, 1-264aa; MGSSHHHHHH SSGLVPRGSH MGSMADSFSL NDALAGSGNP NPQGYPGAWG NQPGAGGYPG AAYPGAYPGQ APPGAYPGQA PPGAYPGQAP PSAYPGPTAP GAYPGPTAPG AYPGSTAPGA FPGQPGAPGA YPSAPGGYPA AGPYGVPAGP LTPYDLPLP GGVMRMLIT IMGTVKPNAN RIVLDFRRGN DVAHFHNPFRF NENRRRIVC NTKQDNNWVGK EERQSAFPFE SGKPFKIQVL VEADHFQVAV NDAHLLQYNH RMKNLREISQ LGISGDITLT SANHAMI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Galectin-3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.