

## Mouse recombinant?Prostatic Acid Phosphatase/ACPP protein

Catalog Number: PROTQ8CE08

### Overview

Product Name	Mouse recombinant?Prostatic Acid Phosphatase/ACPP protein
Description	Mouse recombinant?Prostatic Acid Phosphatase/ACPP protein expressed in Baculovirus with His-Tag. Sequence domain: 32-381aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q8CE08, 32-381aa; KELKFVTLVF RHGDRGPIET FPTDPITSS WPQGFQQLTQ WGMEQHYELG SYIRKRYGRF LNDTYKHDQI  
YIRSTDVRT LMSAMTNLAA LFPPEGISIW NPRLLWQPIP VHTVLSSEDR LLYLPFRDCP RFEELKSETL ESEFLKRLH PYKSFLDTLS SLSGFDDQDL  
FGIWSKVYDP LFCESVHNFT LPSWATEDAM IKLKELSELS LLSLYGIHKQ KEKSRLQGGV LVNEILKNMK LATQPQKYKK LVMYSAHDTT  
VSQLQMALDV YNGVLPPYAS CHMMELYHDK GGHFVEMYR NETQNEPYPL TLPGCTHSCP LEKFAELLDV VISQDWATEC MATSSHQGRN

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant?Prostatic Acid Phosphatase/ACPP protein

For Research Use Only. Not for use in diagnostic procedures.