

## DCUN1D3 (NM\_173475) Human Recombinant Protein

Catalog Number: PROTQ8IWE4

### Overview

Product Name	DCUN1D3 (NM_173475) Human Recombinant Protein
Description	Recombinant protein of human DCN1, defective in cullin neddylation 1, domain containing 3 (S. cerevisiae) (DCUN1D3)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MGQCVTKCKNPSSTLGSKNGDREPSNKSHSRRGAGHREEQVPPCGKPGGDILVNGTKKAEAAATEACQLPTSSGDAGRESKSNAEESLQRLEELFRRYK  
DEREDAILEEGMERFCNDLCVDPTEFRVLLLAWKFQAATMCKFTRKEFFDGCKAISADSIDGICARFPSLLTEAKQEDKFKDLYRFTFQFGLDSEEGQRSL  
HREIAIALWKLVTQNNPPVLDQWLNFLTENPSGIGISRDTWNMFLNFTQVIGPDLSNYSEDEAWPSLFDTFVEWEMERRKREGEGRGALSSGPEGLC  
PEEQT

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



DCUN1D3 (NM\_173475) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.