

Mouse Galectin-4 Recombinant Protein

Catalog Number: PROTQ8K419

Overview

Product Name	Mouse Galectin-4 Recombinant Protein
Description	Mouse Galectin-4 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-326aa. Application(s): SDS-PAGE, Bioactivity.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q8K419, 1-326aa; MAYVPAP GYQPTYNPTL PYKRPIPGGL SVGMSVYIQG MAKENMRRFH VNFAVGQDDG ADVAFHFNPR
FDGWDKVVFN TMQSGQWGKE EKKKSMFQK GKHFEVFMV MPEHYKVVVN GNSFYEGHR LPVQMVTHLQ VDGDLQLSI NFLGGQPAAA
PYPGAMTIPA YPAGSPGYNP PQMNTLPVMT GPPVFNPRVP YVGALQGGLT VRRTHIIKGY VLPTARNFVI NFKVGSSGDI ALHLNPRIGD
SVVRNSFMNG SWGAEERKVA YNPFQPGQFF DLSIRCGMDR FKFVANGQHL FDFSHRFQAF QMVDLEING DITLSYVQI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Galectin-4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.