

Human HSD17B11 Recombinant Protein

Catalog Number: PROTQ8NBQ5

Overview

| | |
|--------------|--|
| Product Name | Human HSD17B11 Recombinant Protein |
| Description | Human HSD17B11 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 20-285aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5mM DTT, 20% glycerol |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q8NBQ5, 20-285aa; MGSSHHHHHH SGLVPRGSH MESFVKLFIP KRRKSVTGEI VLITGAGHGI GRLTAYEFAK LKSKLVLWDI
NKHGLEETAA KCKGLGAKVH TFVVDSCNRE DIYSSAKVKV AEIGDVSILV NNAGVVYTSDFATQDPQIE KTFEVNVLAH FWTTKAFLPA
MTKNNHGHIV TVASAAGHVS VPFLAYCSS KFAAVGFHKT LTDELAALQI TGVKTTCLCP NFVNTGFIKN PSTSLGPTLE PEEVVNRLMH
GILTEQKMIF IPSSIAFLT LERILPERFLAVLKQKI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HSD17B11 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.