

ZNF784 (NM_203374) Human Recombinant Protein

Catalog Number: PROTQ8NCA9

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	ZNF784 (NM_203374) Human Recombinant Protein
Description	Recombinant protein of human zinc finger protein 784 (ZNF784)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MAAARPEAQSRSSPTPESRSQEPLDLVLPDDCRPGTPPSDLIEIQVVKVTDTTLVPEPPEPGSFHCALCPAAFRLVSELLFHEHGHLAGAEGGGQGQGGDP
SRCHVCGHSCPGPASLRAHYSLHTGERPYRCALCPRAFKALAPLLRHQHRHGVEPGTSRRPPDTAAVAEQRPGVAPERAEVVMAAAAAGAAVVGKPF
A
CRFCAKPFRRSSDMRDHERVHTGERPYHCGICGKGFTQSSVLSGHARIHTGERPFRCTLCDRTFNSSNFRKHQRTHFHGPGPGGLGDSGGQLGSSAAE
GSGSGCGVGDPAEEGRGETAKVKVEADQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

ZNF784 (NM_203374) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.