

# Mouse Angiopoietin-like Recombinant Protein 7/ANGPTL7 Recombinant Protein

Catalog Number: PROTQ8R1Q3

## Overview

Product Name	Mouse Angiopoietin-like Recombinant Protein 7/ANGPTL7 Recombinant Protein
Description	Mouse Angiopoietin-like Recombinant Protein 7/ANGPTL7 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 22-337aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

## Concentration

0.25mg/ml (determined by absorbance at 280nm)

## Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

## Purity

> 90% by SDS-PAGE

## Amino Acid Sequence

Mouse, Q8R1Q3, 22-337aa; QKPHKRKTQL KAAGCCEEMR ELKAQVANLS SLLGELSRKQ ESDWVSVMQ VMELESSSKH MESRLSTAES  
KYSEMNNQID IMQLQAAQTV TQTSADAIYD CSSLYQKNYR ISGVYKLPPD EFLGSPELEV FCDMETSGGG WTIIQRRKSG LVSFYQDWRQ  
YKQGFSGIRG DFWLGNEHIH RLTRQPSRLR VELEDWEGNA RYAEYSYFAL GNELNSYRLF LGNYSGNVGK DALLYHNNTV FSTKDKDNDN  
CLDKCAQLRK GGYWYNCCTD SNLNGVYYRL GEHRKHMDGI SWYGWHGANY SLKRVEMKIR PEAFTP ?/description\_after\_attributes?

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Angiopoietin-like Recombinant Protein 7/ANGPTL7 Recombinant Protein