

## Human TIM-3 Recombinant Protein

Catalog Number: PROTQ8TDQ0

### Overview

Product Name	Human TIM-3 Recombinant Protein
Description	Human TIM-3 Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 22-202aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q8TDQ0, 22-202aa; ADLSEVEYRA EVGQNAYLPC FYTPAAPGNL VPVCWVGKGAC PVFECGNVVL RTDERDVNYW TSTRYWLNGDF  
RKGDVSLTIE NVTLADSGIY CCRIQIPGIM NDEKFNLKLK IKPAKVTPAP TRQRDFTAAF PRMLTTRGHG PAETQTLGSL PDINLTQIST LANELRDSRL  
ANDLRDSGAT IRIGLEPKSC DKTHCPCPCP APELLGGPSV FLFPPKPKDT LMISRTPEVT CVVVDVSHED PEVKFNWYVD GVEVHNAKTK  
PREEQYNSTY RVVSVLTVLH QDWLNGKEYK CKVSNKALPA PIEKTISKAK GQPREPQVYV LPPSRDELTK NQVSLTCLVK GFYPSDIAVE  
WESNGQPENN YKTTTPVLDS DGSFFLYSKL TVDKSRWQQG NVFSCSVME ALHNHYTQKS LSLSPGKHHH HHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human TIM-3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.