

## PP4R1 (PPP4R1) (NM\_005134) Human Recombinant Protein

Catalog Number: PROTQ8TF05

### Overview

Product Name	PP4R1 (PPP4R1) (NM_005134) Human Recombinant Protein
Description	Recombinant protein of human protein phosphatase 4, regulatory subunit 1 (PPP4R1), transcript variant 2
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Concentration

>50 µg/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MADLSLLQEDLQEDADGSLDFVSQDEMLTPLGRDLKYAASENIFNRQMVARSLDLTREVCDDEDCIAVLERISRLADDSEPTVRAELMEQVPHIALFC  
QENRPSIPYAFSKFLLPIVVRYLADQNNQVRKTSQAALLALLEQELIERFDVETKVCPLIELTAPDSNDDVKTEAVAIMCKMAPMVGKDITERLILPRFCE  
MCCDCRMFHVVRKVCANFGDICSVVGQQATEEMLLPRFFQLCSDNVWGVKACAECFMAVSCATCQEIRRTKLSALFINLISDPSRWVRQAAFQSLG  
PFISTFANPSSSGQYFKEESKSSEEMSVENKNRTRDQEAPEDVQVRPEDTPSDLSVSNSSVILENTMEDHAAEASGKPLGEISVPLDSSLLCTLSSESHQEA  
ASNENDKKPGNYKSMRLPEVGTTSQDSALLDQELYNFHFWRTPLEIDLIDIELEQNSGGKPSPEGPEEESEGPVPPSPNITMATRKELEEMIENLEPHID  
DPDVKAQVEVLSAALRASSLDAHEETISIEKRSDLQDELINELPNCKINQEDSVPLISDAVENMDSTLHYIHSDSDLSNNSFSFPDEERRTKVQDVVPQAL  
LDQYLSMTDPSRAQTVDEIAKHCAYSPLGVALTLGRQNWHLCLRETYETLASDMQWKVRRTLAFSIHELAVILGDQLTAADLVPIFNGFLKDLDEVRIGV  
LKHLHDFLKLHIDKRREYLYQLQEFVLDNSRNWRFRAELAEQLLILLELYSPRDVYDYLRPIALNLCADKVSSVRWISYKLVSEMVKKLHAATPPTFGVD  
LINELVENFGRCPKWSGRQAFVFCQTVIEDDCLPMDQFVHLMPHLLTLANDRVPNVVLLAKTLRQTLLKEDYFLASASCHQEAVEQTIMALQMDR  
DSDVKYFASIH PASTKISEDAMSTASSTY

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



PP4R1 (PPP4R1) (NM\_005134) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.