

Mouse recombinant TIM-3/HAVCR2 protein

Catalog Number: PROTQ8VIM0

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant TIM-3/HAVCR2 protein
Description	Mouse recombinant TIM-3/HAVCR2 protein expressed in HEK293 with hIgG-His-Tag. Sequence domain: 20-193aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Mouse, Q8VIM0, 20-193aa; RSLNAYVFE VGKNAYLPCS YTLSTPGALV PMCWVGKGFCP WSQCTNELLR TDERNVTYQK SSRYQLKGDL NKGDVSLIIK NVTLDHGTY CCRIQFPLM NDKKLELKLD IKAQKVTAPQ TAHGDSSTAS PRTLTERNG SETQTLVTLH NNGTKISTW ADEIKDSGET IRTA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant TIM-3/HAVCR2 protein](#)

For Research Use Only. Not for use in diagnostic procedures.