

GTF3C2 (NM_001521) Human Recombinant Protein

Catalog Number: PROTQ8WUA4

Overview

Product Name	GTF3C2 (NM_001521) Human Recombinant Protein
Description	Recombinant protein of human general transcription factor IIIC, polypeptide 2, beta 110kDa (GTF3C2), transcript variant 1
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MDTCGVGYVALGEAGPVGNM TVVDS PGQEV LNQLDVKTSSEMTSAEASVEMSLPTPLPGFEDSPDQRRLLPPEQESLSRLEQPDLSSSEMSKVKPRASKP
GRKRGRTRKGPKRPPQPNPPSAPLVPGLLDQSNPLSTPMPKKRGRKSKAE LLLLKSKDLDRPESQSPKRPPEDFETPSGERPRRRRAAQVALLYLQELA
EELSTALPAPVSCPEGPKVSSPTKPKKIRQPAACPGGEEVDGAPRDEDFLQVEAEDVEESEGPSESSSEPEPVVPRSTPRGSTSGKQKPHCRGMAPNGL
PDHIMAPVWKCLHLTKDFREQKHSYWEFAEWIPLAWKWHLLELEAAPYLPQEEKSPLFSVQREGLPEDGTLYRINRFSSITAHPERWDVSFFTGGPLV
ALDWCPVPEGAGASQYVALFSSPDMNETHPLSQLHSGPGLLQLWGLGTLQQESC PGNRAHFVYGIACDNGCIWDLKFCPSGAWELPGTTPRKAPLLPR
LGLLALACSDGKVLFLSLPHPEALLAQPPDAVKPAIYKVCVATLQV GSMQATDPSECGQCLSLAWMPTRPHQHLAAGYYNGMVVFWNLPTNSPLQ
RIRLSDGSLKLYPFQCFLAHDQAVRTLQWCKANSHFLVSAGSDRKIKFWDLRRPYEPINSIKRFLSTELAWLLPYNGVTVAQDNCYASYGLCGIHYIDAG
YLGFKAYFTAPRKGTVWLSGSDWLGTIAAGDISGELIAAILPDMALNPINVKRPVERRFPIYKADLIPYQDSPEGPDHSSASSGVPNPPKARTYTETVNH
HYLLFQD TDLGSFHDLRREPLMRMQEGEGHSQLCLDRLQLEAIHKVRFSPNLDSYGWLVS GGQSLVRIHFVRGLASPLGHRMQLESRAHFNAMEFQP
SSPTRRPGFSPTSHRLLPTP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



GTF3C2 (NM_001521) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.