

Human RAB2B Recombinant Protein

Catalog Number: PROTQ8WUD1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human RAB2B Recombinant Protein
Description	Human RAB2B Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-216aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT, 0.15M NaCl.

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q8WUD1, 1-216aa; MTYAYL FKYIIIGDTG VGKSCLLQF TDKRFQPVHD LTIGVEFGAR MVNIDGKQIK LQIWDTAGQE SFRSITRSYY
RGAAGALLVY DITRRETFNH LSWLEDARQ HSSSNMVI MLIGNKSDLES RDVKREEGEA FAREHGLIFM ETSAKTACNV EEFINTAKE
IYRKIQGGLF DVHNEANGIK IGPQQSISTS VGPSASQRNS RDIGSNSGCC

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human RAB2B Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.