

Stabilin 2 (STAB2) (NM_017564) Human Recombinant Protein

Catalog Number: PROTQ8WWQ8

Overview

| | |
|--------------|--|
| Product Name | Stabilin 2 (STAB2) (NM_017564) Human Recombinant Protein |
| Description | Recombinant protein of human stabilin 2 (STAB2) |
| Source | HEK293T |
| Formulation | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol |

Storage

Store at -80°C.

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Stabilin 2 (STAB2) (NM_017564) Human Recombinant Protein