

## Mouse Prss28 Recombinant Protein

Catalog Number: PROTQ924N9

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                                                                                               |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Mouse Prss28 Recombinant Protein                                                                                                                                                              |
| Description  | Mouse Prss28 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 27-274aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size         | Starting from 10ug                                                                                                                                                                            |
| Tag          | His-Tag                                                                                                                                                                                       |
| Form         | Liquid                                                                                                                                                                                        |
| Source       | Baculovirus                                                                                                                                                                                   |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol                                                                                                                                    |

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q924N9, 27-274aa; KPVGIVGGQC TPPGKWPWQV SLRMYSYEVN SWVHICGCSI IHPQWILTAA HCIQSQDADP AVYRVQVGEV  
YLYKEQELLN ISRIIHPDY NDVSKRFDLA LMQLTALLVT STNVSPVSLP KDSSTFDSTD QCWLVGWGNL LQRVPLQPPY QLHEVKIPIQ  
DNKSKRAYR KKSSDEHKAV AIFDDMLCAG TSGRGPCFGD SGGPLVCWKS NKWIQVGVVS KGIDCSNNLP SIFSRVQSSL AWIHQHIQLE  
HHHHHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Mouse Prss28 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.