

SEC14L1 (NM_003003) Human Recombinant Protein

Catalog Number: PROTQ92503

Overview

| | |
|--------------|--|
| Product Name | SEC14L1 (NM_003003) Human Recombinant Protein |
| Description | Recombinant protein of human SEC14-like 1 (<i>S. cerevisiae</i>) (SEC14L1), transcript variant 1 |
| Size | 20 µg |
| Tag | C-Myc/DDK |
| Form | Frozen Solution in PBS Buffer |
| Source | HEK293T |
| Formulation | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol |

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MVQKYQSPVRVYKYPFELIMAAAYERRFPTCLIPMFVGS DTVNEFKSE DGAIHVIERRCKLDVDAPRL LKKIAGVDYVYFVQKNSLNSRRTLHIEAYNET
FSNRVIINEHCCYTVHPENEDWTCFEQSASLDIKSFFGFESTVEKIAMKQYTSNIKKGKEIIEYLRQLEEEGITFVPRWSPPSITPSETSSSSSKKQAASMA
VVIPEAALKEGLSGDALSSPSAPEPVVGT PDDKLDADYIKRYLGDLTPLQESCLIRLRQWLQETHKGIKPKDEHILRFLRARDFNIDKAREIMCQSLTWRK
QHQVDYILETWTPPQVLQDY YAGGWHHHDKGRPLYVLR LQGMDTKGLVLRALGEEALLRYVLSINEEGLRRCEENTKVFGRPISSWTCLVDLEGLNM
RHLWRPGVKALLRIIEVVEANYPETLGRLLILRAPRVFPVLWTLVSPFIDDNTRRKFLIYAGNDYQGPGLLDYIDKEIIPDFLSGECMCEVPEGGLVPKSLY
RTAEELNEDLKLWTETIYQSASVFKGAPHEILIQIVDASSVITWDFDVCKGDIVFNIYHSKRSPQPPKKDSLGAHSITSPGGNNVQLIDKVVQLGRDYSM
VESPLICKEGESVQGS HVTRWPGFYILQWKFHSM PACAASSLPRVDDVLASLQVSSHKCKVMYYTEVIGSEDFRGSMTSLESSHSGFSQLSAATTSSSQS
HSSSMISR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



SEC14L1 (NM_003003) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.