

## NSMAF (NM\_003580) Human Recombinant Protein

Catalog Number: PROTQ92636

### Overview

Product Name	NSMAF (NM_003580) Human Recombinant Protein
Description	Recombinant protein of human neutral sphingomyelinase (N-SMase) activation associated factor (NSMAF), transcript variant 1
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MAFIRKKQQEQQLQLYSKERFSLLLLNL E E Y F E Q H R A N H I L H K G S H H E R K I R G S L K I C S K S V I F E P D S I S Q P I I K I P L R D C I K I G K H G E N G A N R H F T K A K S G G  
I S L I F S Q V Y F I K E H N V V A P Y K I E R G K M E Y V F E L D V P G K V E D V V E T L L Q L H R A S C L D K L G D Q T A M I T A I L Q S R L A R T S F D K N R F Q N I S E K L H M E C K A E M V T P  
L V T N P G H V C I T D T N L Y F Q P L N G Y P K P V V Q I T L Q D V R R I Y K R R H G L M P L G L E V F C T E D D L C S D I Y L K F Y E P Q D R D D L Y F Y I A T Y L E H H V A E H T A E S Y M L Q W  
Q R G H L S N Y Q Y L L H L N N L A D R S C N D L S Q Y P V F P W I I H D Y S S E L D L S N P G T F R D L S K P V G A L N K E R L E R L L T R Y Q E M P E P K F M Y G S H Y S S P G Y V L F Y L V R I A  
P E Y M L C L Q N G R F D N A D R M F N S I A E T W K N C L D G A T D F K E L I P E F Y G D D V S F L V N S L K L D L G K R Q G G Q M V D D V E L P P W A S S P E D F L Q K S K D A L E S N Y V S  
E H L H E W I D L I F G Y K Q K G S D A V G A H N V F H P L T Y E G G V D L N S I Q D P D E K V A M L T Q I L E F G Q T P K Q L F V T P H P R R I T P K F K S L S Q T S S Y N A S M A D S P G E E S F E  
D L T E E S K T L A W N N I T K L Q L H E H Y K I H K E A V T G I T V S R N G S S V F T T S Q D S T L K M F S K E S K M L Q R S I S F S N M A L S S C L L L P G D A T V I T S S W D N N V Y F Y S I A F G R  
R Q D T L M G H D D A V S K I C W H D N R L Y S A S W D S T V K V W S G V P A E M P G T K R H H F D L L A E L H D V S V D T I S L N A A S T L L V S G T K E G T V N I W D L T T A T L M H Q I P  
C H S G I V C D T A F S P D S R H V L S T G T D G C L N V I D V Q T G M L I S S M T S D E P Q R C F V W D G N S V L S G S Q S G E L L V W D L L G A K I S E R I Q G H T G A V T C I W M N E Q C S S I I  
T G G E D R Q I I F W K L Q Y

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



NSMAF (NM\_003580) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.