

Human ART4 Recombinant Protein

Catalog Number: PROTQ93070

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human ART4 Recombinant Protein
Description	Human ART4 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 47-285aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q93070, 47-285aa; SEVAIKI DFDFAFGSFD DQYQGCSKQV MEKLTQGDYF TKDIEAQKNY FRMWQKAHLA WLNQGGKVLPPQ NMTTTHAVAI LFYTLNSNVH SDFTRAMASV ARTPQQYERS FHFYLYHYLL TSAIQLLRKD SIMENGTLCY EVHYRTKDVH FNAYTGATIR FGQFLSTSL KEEAQEFGNQ TLFTIFTCLG APVQYFSLKK EVLIPPYELF KVINMSYHPR GDWLQLRSTG NLSTYNCQLL KA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human ART4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.