

## Canine recombinant IL-5 protein

Catalog Number: PROTQ95J76

### Overview

|              |   |
|--------------|---|
| Product Name | Canine recombinant IL-5 protein   |
| Description  | Canine recombinant IL-5 protein expressed in Baculovirus with His-Tag. Sequence domain: 22-134aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size         | Starting from 10ug  |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | Baculovirus   |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol  |

### Concentration

0.25mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Canine, Q95J76, 22-134aa; VENPMNRLVA ETLTLLSTHR TWLIGDGNLM IPTPENKNHQ LCIKEVFQGI DTLKNQTAHG EAVDKLFQNL SLIKEHIERQ KKRCAGERWR VTKFLDYLVQV FLGVINTEWT PES

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Canine recombinant IL-5 protein

For Research Use Only. Not for use in diagnostic procedures.