

## Human UBE2E3 Recombinant Protein

Catalog Number: PROTQ969T4

### Overview

Product Name	Human UBE2E3 Recombinant Protein
Description	Human UBE2E3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-207aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol, 1mM DTT

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 80% by SDS-PAGE

### Amino Acid Sequence

Human, Q969T4, 1-207aa; MGSSHHHHHH SSGLVPRGSH MGSMSSDRQR SDDESPSTSS GSSDADQRDP AAPEPEEQEE RKPSATQQKK  
NTKLSKTTA KLSTSAKRIQ KELAEITLDP PPNC SAGPKG DNIYEW RSTI LGPPGSVYEG GVFFLDITFS SDYPFKPPKV TFRTRIYHCN INSQGVICLD  
ILKDNWSPAL TISKVLLSIC SLLTDCNPAD PLVGS IATQY LTNRAEHDRI ARQWTKRYAT

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human UBE2E3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.