

## Human MAGOHB Recombinant Protein

Catalog Number: PROTQ96A72

### Overview

|              |   |
|--------------|---|
| Product Name | Human MAGOHB Recombinant Protein  |
| Description  | Human MAGOHB Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-148aa. Application(s): SDS-PAGE. |
| Size         | Starting from 10ug  |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | E. coli   |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT 30% glycerol, 0.15M NaCl   |

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q96A72, 1-148aa; MGSSHHHHHH SGLVPRGSH MGSMASDF YLRYVGHKG KFGHEFLEFE FRPDGKLRYA NNSNYKNDVM IRKEAYVHKS VMEELKRIID DSEITKEDDA LWPPPDRVGR QELEIVIGDE HISFTTSKIG SLIDVNQSKD PEGLRVFYLL VQDLKCLVFS LIGLHFKIKP I

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human MAGOHB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.