

Human ACSF2 Recombinant Protein

Catalog Number: PROTQ96CM8

Overview

Product Name	Human ACSF2 Recombinant Protein
Description	Human ACSF2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 42-615aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q96CM8, 42-615aa; LSSREVD RMVSTPIGGL SYVQGCTKKH LNSKTVGQCL ETTAQRVPER EALVVLHEDV RLTFACLKEE VDKAASGLLS IGLCKGDRLG MWGPNSYAWW LMQLATAQAG IILVSVNPAY QAMELEYVLK KVGCKALVFP KQFKTQQYYN VLKQICPEVE NAQPGALKSQ RLPDLTTVIS VDAPLPGTLL LDEVVAAGST RQHLDQLQYN QQFLSCHDPI NIQFTSGTTG SPKGATLSHY NIVNNSNILG ERLKLHEKTP EQLRMILPNP LYHCLGSVAG TMMCLMYGAT LILASPIFNG KKALEAISRE RGTFLYGTPT MFVDILNQP D FSSYDISTMC GGVIAGSPAP PELIRAIINK INMKDLVVAY GTTENSPTVF AHFPEDTVEQ KAESVGRIMP HTEARIMNME AGTLAKLNTP GELCIRGYCV MLGYWGEPQK TEEAVDQDKW YWTGDVATMN EQGFCKIVGR SKDMIIRGGE NIYPAELEDF FHHPKVQEV QVVGKDDRM GEEICACIRL KDGEETTVEE IKAFCKGKIS HFKIPKYIVF VTNYPLTISG KIQKFKLREQ MERHLNL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ACSF2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.