

R3HDM4 (NM_138774) Human Recombinant Protein

Catalog Number: PROTQ96D70

Overview

| | |
|--------------|---|
| Product Name | R3HDM4 (NM_138774) Human Recombinant Protein |
| Description | Recombinant protein of human chromosome 19 open reading frame 22 (C19orf22) |
| Size | 20 µg |
| Tag | C-Myc/DDK |
| Form | Frozen Solution in PBS Buffer |
| Source | HEK293T |
| Formulation | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol |

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MVALENPECGPEAAEGTPGGRRLLPLPSCLPALASSQVKRLSASRRKQHFINQAVRNSDLVPAKAKGRKSLQRLENTQYLLTLETDGGLPGLEDGDLAPP
ASPGIFAEACNNATYVEVWDFMNRSGEEQERVLRYLEDEGRSKARRRGPRGEDRRREDPAYTPRECFQRISRRLRAVLKRSRIPMETLETWEERLLRF
FSVSPQAVYTAMLDNSFERLLLHAVCQYMDLISASADLEGKRQMKVSNRHLDLPPGLLLSAYLEQHS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



R3HDM4 (NM_138774) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.