

## Human Musashi-2 Recombinant Protein

Catalog Number: PROTQ96DH6

### Overview

Product Name	Human Musashi-2 Recombinant Protein
Description	Human Musashi-2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-328aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q96DH6, 1-328aa; MEANGSQ GTSGSANDSQ HDPGKMFIGN LSWQTSPDSL RDYFSKFGEI RECMVMRDPT TKRSRFGFV  
TFADPASVDK VLGQPHHELD SKTIDPKVAF PRRAQPKMVT RTKKIFVGGI SANTVVEDVK QYFEQFGKVE DAMLMFDKTT NRHRGFGFVT  
FENEDVVEKV CEIHFHEINN KMVECKKAQP KEVMFPPGTR GRARGLPYTM DAFMLGMGML GYPNFVATYG RGYPGFAPSY GYQFPGFPAA  
AYGPVAAA AV AAARGSGSNP ARPGGFPGAN SPGPVADLYG PASQDSGVGN YISAASPQPG SGFGHGIAGP LIATAFTNGY H

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Musashi-2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.