

Human Dysbindin Recombinant Protein

Catalog Number: PROTQ96EV8

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human Dysbindin Recombinant Protein
Description	Human Dysbindin Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-270aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 0.5mM DTT, 20% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q96EV8, 1-270aa; MLSA HWEKKKTSLV ELQEQLQQLP ALIADLESMT ANLTHLEASF EEVENNLLHL EDLCGQCELE RCKHMQSQQ
LENYKKNKRKE LETFKAELDA EHAQKVLEME HTQQMKLKER QKFFEEAFQQ DMEQYLSTGY LQIAERREPI GSMSSMEVNV DMLEQMDLMD
ISDQEALDVF LNSGGEENTV LSPALGPRESS TCQNEITLQV PNPSELRAKP PSSSSTCTDS ATRDISEGGE SPVVQSDEEE VQVDTALATS
HTDREATPDG GEDSDS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human Dysbindin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.