

## CGI 62 (ZC2HC1A) (NM\_016010) Human Recombinant Protein

Catalog Number: PROTQ96GY0

### Overview

Product Name	CGI 62 (ZC2HC1A) (NM_016010) Human Recombinant Protein
Description	Recombinant protein of human family with sequence similarity 164, member A (FAM164A)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MEGLEENGGVVQVGELLPCKICGRFFPVALKKHGPICQKTATKKRKTFDSSRQRAEGTDIPTVKPLKPRPEPPKPSNWRRKH EEFIATIRAAKGLDQA  
LKEGGKLP P P P P P S Y D P D Y I Q C P Y C Q R R F N E N A A D R H I N F C K E Q A A R I S N K G K F S A D T K G K P T S R T Q V Y K P P A L K K S N S P G T A S S G S S R L P Q P S G A G K T V  
V G V P S G K V S S S S S L G N K L Q T L S P S H K G I A A P H A G A N V K P R N S T P P S L A R N P A P G V L T N K R K T Y E S Y I A R P D G C A S S L N G G N I K G I E G H S P G N L P K F C H  
E C G T K Y P V E W A K F C C E C G I R R M I L

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



CGI 62 (ZC2HC1A) (NM\_016010) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.