

## LRIG1 (NM\_015541) Human Recombinant Protein

Catalog Number: PROTQ96JA1

### Overview

|              |   |
|--------------|---|
| Product Name | LRIG1 (NM_015541) Human Recombinant Protein   |
| Description  | Recombinant protein of human leucine-rich repeats and immunoglobulin-like domains 1 (LRIG1) |
| Size         | 20 µg   |
| Tag          | C-Myc/DDK   |
| Form         | Frozen Solution in PBS Buffer   |
| Source       | HEK293T   |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol  |

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MARPVRGGLGAPRRSPCLLLLWLLLLRLEPVTAAAGPRAPCAAACCTCAGDSLDCGGRGLAALPGDLPWSTRSLNLSYNKLEIDPAGFEDLPNLQEVYLN  
NNELTAVPSLGAASSHVVSFLQHNKIRSVESQLKAYLSLEVLDLSLNNITEVRNTCFPHGPPIKELNLAGNRIGTLELGAFDGLSRSLTLRLSKNRITQLP  
VRAFKLPRLTQLDLNRRNRIEGLTFQGLNSLEVLKLRNNISKLTDFGAFWGLSKMHVLHLEYNSLVEVNSGSLYGLTALHQLHLSNNSIARIHRKQWSF  
CQKLHELVLFSFNNLRLDEESLAELSSLSVLRSHNSISHIAEGAFKGLRSLRVLDLDHNEISGTIEDTSGAFSGLDSLKLTFLGNIKIKSVAKRAFSGLEGLEH  
LNLGGNAIRSVQFADFVKMKNLKEHISSDSFLCDCQLKWLPPWLI GRMLQAFVTATCAHPESLKGQSIFSVPPESFVCDLKPQIITQPETTMAMVVGK  
DIRFTCSAASSSSPMTFAWKDNEVLTNADMENFVHVHAQDGEVMEYTTILHLRQVTFGHEGRYQCIVITNHFGSTYSHKARLTVNVLPSFTKTPHDIT  
IRTTT MARLECAATGHPNPQIAWQKDGDTDFPAARERRMHVMPDDDDVFFITDVKIDDAGVYSCTAQNSAGSISANATLTVLETPSLVVPLEDRVVS VG  
ETVALQCKATGNPPPRITWFKGDRPLSLTERHHLTPDNQLLVVQNVVAEDAGRYTCEMSNTLGTERRAHSQSVLPAAGCRKDGTTVGIFTIAVSSIVLT  
SLVWVCIYQTRKKSEEYSVTNTDETVPDPVPSYLSSQGTLSDRQETVVRTEGGPQANGHIESNGVCPRDASHFPEPDTHSVACRQPKLCAGSAYHKE  
PWKAMEKAEGTPGPHKMEHGGRVVCSDCNTEVDCYSRQAFHPQPVS RDSAQPSAPNGPEPGGSDQEHSPHHQCSRTAAGSCPECQGSLYPSNHD  
RMLTAVKKKPMASLDGKGDSSWTLARLYHPDSTELQPASLTSGSPERA EAQYLLVSNHLPKACDASPESTPLTGQLPGKQRVPLLLAPKS

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



LRIG1 (NM\_015541) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.