

## Human ZNF514 Recombinant Protein

Catalog Number: PROTQ96K75

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                                                                                                                                                                              |
|--------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Human ZNF514 Recombinant Protein                                                                                                                                                                                                                                             |
| Description  | Human ZNF514 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-400aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays. |
| Size         | Starting from 20ug                                                                                                                                                                                                                                                           |
| Tag          | His-Tag                                                                                                                                                                                                                                                                      |
| Form         | Liquid                                                                                                                                                                                                                                                                       |
| Source       | E. coli                                                                                                                                                                                                                                                                      |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol                                                                                                                                                                                                             |

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q96K75, 1-400aa; MGSSHHHHHH SSGLVPRGSH MGSMTFEDVA VEF SQWEWGQ LNPAQKDLYR EVMLNFRNL AILG LLVSKP YVICQLEEG EPFMVEREIS TGAHSDWKRR SKSKESMP SW GISKEELFQV VSVEKHIQDV LQFSKLKAAC GCDGQLEM QQ IKQERHLKQM STIHK SATTL SRDYKWNGFG RSLGLRSVLV NQHSILMGE G SYKCDTEFRQ TLGGNNSQRT HPEKKSCKCN ECGKSFHFQS ELRRHQ RCHT GEKPYECSDC GRAFGHISL IKHQ RTH TGE KPYECSECGR AFSQSSSLVL HYRFHTGEKP YKCNECGRAF GHTSSLIKHQ RTH TGEKPYE

CRECGRTFSQ SSSLIVHYRF HTGEKPYKCN KCGRAFSQSS SLTQHYRFHT GEKPYKCNEC GRAFAHTASL IKHQRSHAGK KTL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### [Human ZNF514 Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.