

Human CNDP2/CPGL Recombinant Protein

Catalog Number: PROTQ96KP4

Overview

Product Name	Human CNDP2/CPGL Recombinant Protein
Description	Human CNDP2/CPGL Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-475aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q96KP4, 1-475aa; MGSSHHHHHH SGLVPRGSH MGSMAALTTL FKYIDENQDR YIKKLAKWVA IQSVSAWPEK RGEIRRMMEV
AAADVQLGG SVELVDIGKQ KLPDGSEIPL PPILLGRLGS DPQKKTVCYI GHLDVQPAAL EDGWDSEPF LVERDGKLYG RGSTDDKGPV
AGWINALEAY QKTGQEIPVN VRFCLEGMEE SGSEGLDELI FARKDTFFKD VDYVCISDNY WLGKKKPCIT YGLRGICYFF IEVECSNKDL
HSGVYGGSVH EAMTDLILLM GSLVDKRGNI LIPGINEAVA AVTEEEHKLY DDIDFDIEEF AKDVGAQILL HSHKKDILMH RWRYPSSLH
GIEGAFSGSG AKTVIPRKVV GKFSIRLVPN MTPEVVGEQV TSYLTKKFAE LRSNPEFKVY MGHGGKPVWS DFSHPHYLAG RRAMKTVFGV
EPDLTREGGS IPVTTLTFQEA TGKNVMLLPV GSADDGAHSQ NEKLNRYNYI EGTKMLAAYL YEVSQDKD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CNDP2/CPGL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.