

## Human Siglec-10 Recombinant Protein

Catalog Number: PROTQ96LC7

### Overview

Product Name	Human Siglec-10 Recombinant Protein
Description	Human Siglec-10 Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 17-455aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q96LC7, 17-455aa; MDGRFWIRVQ ESVMVPEGLC ISVPCSFSSYP RQDWTGSTPA YGYWFKAVTE TTKGAPVATN HQSREVEMST RGRFQLTGDP AKGNCSLVIR DAQMQDESQY FFRVERGSYV RYNFMNDGFF LKVTALTQKP DVYIPETLEP GQPVTVICVF NWAFFEECPPP SFSWTGAALS SQGTKPTTSH FSVLSFTPRP QDHNTDLTCH VDFSRKGVSA QRTVRLRVAY APRDLVISIS RDNTPALEPQ PQGNVPYLEA QKGQFLRLLC AADSQPPATL SWVLQNRVLS SSWPWGPRPL GLELPGVKAG DSGRYTCRAE NRLGSQQRAL DLSVQYPPEP LRVMVSQLNR TVLENLNGNT SLPVLEGQSL CLVCVTHSSP PARLSWTQRG QVLSPSQPSD PGVLELPRVQ VEHEGEFTCH ARHPLGSQHV SLSLSVHYKK GLISTAFSN

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Siglec-10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.