

Human RDH12 Recombinant Protein

Catalog Number: PROTQ96NR8

Overview

Product Name	Human RDH12 Recombinant Protein
Description	Human RDH12 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 39-316aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 40% glycerol, 0.2M NaCl

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q96NR8, 39-316aa; GKVVV ITGANTGIGK ETARELASRG ARVYIACRDV LKGESAASEI RVDTKNSQVL VRKLDLSDTK SIRAF AEGFL AEEKQLHILI NNAGVMMCPY SKTADGFETH LGVNH LGHFL LTYLLERLK VSAPARVVNV SVAHHIGKI PFHDLQSEKR YSRGFAYCHS KLANVLF TRE LAKRLQGTGV TTYAVHPGVV RSELVRHSSL LCLLWRLFSP FVKTAREGAQ TSLH CALAEG LEPLSGKYFS DCKRTWVSPR ARNNKTAERL WNVSCCELLGI RWE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RDH12 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.