

Human TRIB3 Recombinant Protein

Catalog Number: PROTQ96RU7

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human TRIB3 Recombinant Protein
Description	Human TRIB3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-358aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q96RU7, 1-358aa; MGSSHHHHHH SGLVPRGSH MGSMRATPLA APAGSLSRKK RLEDDNLDLTPVQKRARS GPQRLPPCL
LPLSPPTAPD RATAVATASR LGPYVLEPE EGGRAYQALH CPTGTEYTCK VYPVQEALAV LEPYARLPPH KHVARPTEVL AGTQLLYAFF
TRTHGDMHSL VRSRHRIPPEA AVLFRQMA TALAHCHQHG LVLRDLKLCR VFVADRERKK LVLENLEDSC VLTGPDDSLW DKHACPAYVG
PEILSSRASY SGKAADVWSL GVALFTMLAG HYPFQDSEPV LLFGKIRRGALYPAGLSAP ARCLVRCLLR REPAERLTAT GILLHPWLRQ
DPMPLAPTRS HLWEAAQVVP DGLGLDEARE EEGDREVVLYG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human TRIB3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.