

## Human AMELX Recombinant Protein

Catalog Number: PROTQ99217

### Overview

Product Name	Human AMELX Recombinant Protein
Description	Human AMELX Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 17-191aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.5) containing 0.2M NaCl, 30% glycerol, 1mM DTT

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q99217, 17-191aa; MPLPPHP GHPGYINFSY EVLTPLKWYQ SIRPPYPSYG YEPMGGWLHH QIIPVLSQQH PPTH TLQPHH HIPVVPAQQP  
VIPQQPMMPV PGQHSMTPIQ HHQPNLPPPA QQPYPQPVQ PQPHQPMQPQ PPVHPMQPLP PQLPPLPMFP MQPLPMLPD LTLEAWPSTD  
KTKREEVD

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human AMELX Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.