

VAT1 (NM_006373) Human Recombinant Protein

Catalog Number: PROTQ99536

Overview

Product Name	VAT1 (NM_006373) Human Recombinant Protein
Description	Recombinant protein of human vesicle amine transport protein 1 homolog (T. californica) (VAT1)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 µg/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MSDEREVAEAATGEDASSPPKTEAASDPQHAAASEGAAAAAASPPLLRCLVLTGFGGYDKVKLQSRPAAPPAPGPGQLTLRLRACGLNFADLMARQG
LYDRLPPLPVTPGMEGAGVVIIVGEGVSDRKAGDRVMVLNRSGMWQE E V T V P S V Q T F L I P E A M T F E E A A A L L V N Y I T A Y M V L F D F G N L Q P G H S V L V H
MAAGGVGMAAVQLCRTVENVTVFGTASASKHEALKENG V T H P I D Y H T T D Y V D E I K K I S P K G V D I V M D P L G G S D T A K G Y N L L K P M G K V V T Y G M A N L L T
GPKRNLMA L A R T W W N Q F S V T A L Q L L Q A N R A V C G F H L G Y L D G E V E L V S G V V A R L L A L Y N Q G H I K P H I D S V W P F E K V A D A M K Q M Q E K K N V G K V L L V P
GPEKEN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



VAT1 (NM_006373) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.