

Human Monoglyceride Lipase/MGLL Recombinant Protein

Catalog Number: PROTQ99685

Overview

Product Name	Human Monoglyceride Lipase/MGLL Recombinant Protein
Description	Human Monoglyceride Lipase/MGLL Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-313aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q99685, 1-313aa; METGPEDPSS MPEESSPRRT PQSIPYQDLP HLVNADGQYL FCRYWKPTGT PKALIFVSHG AGEHSGRYEE
LARMLMGLDL LVFAHDHVGH GQSEGERMVV SDFHVFVRDV LQHVDMSMQKD YPGLPVFLLG HSMGGAIAIL TAAERPGHFA GMVLISPLVL
ANPESATTFK VLAACKVNLV LPNLSLGPID SSVLSRNKTE VDIYNSDPLI CRAGLKVCFG IQLLNAVSRV ERALPKLTVP FLLQGSADR LCDSKGAYLL
MELAKSQDKT LKIYEGAYHV LHKELPEVTN SVFHEINMWV SQRTATAGTA SPP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Monoglyceride Lipase/MGLL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.