

Human CCL19/MIP-3 beta Recombinant Protein

Catalog Number: PROTQ99731

Overview

Product Name	Human CCL19/MIP-3 beta Recombinant Protein
Description	Human CCL19/MIP-3 beta Recombinant Protein expressed in E. coli with T7-tag. Sequence domain: 22-98aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	T7-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q99731, 22-98aa; GTND AEDCCLSVTQ KPIPGYIVRN FHYLLIKDGC RVPVAVFTTL RGRQLCAPPD QPWVERIIQR LQRSAKMKR RSS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CCL19/MIP-3 beta Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.