

Human TSG101 Recombinant Protein

Catalog Number: PROTQ99816

Overview

Product Name	Human TSG101 Recombinant Protein
Description	Human TSG101 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-145aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q99816, 1-145aa; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMVAESQLKKMVSK YKYRDLTVRE TVNVITLYKD
LKPVLDSYVF NDGSSRELMN LTGTIPVPYR GNTYNIPICL WLLDTYPYNP PICFVKPTSS MTIKTGKHVD ANGKIYLPYL HEWKHPQSDL
LGLIQVMIVV FGDEPPVFSR P

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human TSG101 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.