

Human BAG1 Recombinant Protein

Catalog Number: PROTQ99933

Overview

Product Name	Human BAG1 Recombinant Protein
Description	Human BAG1 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-230aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 100mM NaCl, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q99933, 1-230aa; MNRSQEVTRD EESTRSEEVT REEMAAAGLT VTVTHSNEKH DLHVTSQQGS SEPVVQDLAQ VVEEVIGVPQ SFQKLIFKGGK SLKEMETPLS ALGIQDGCRV MLIGKKNPQ EEVELKCLKH LEKSVEKIAD QLEELNKELT GIQQGFLPKD LQAEALCKLD RRVKATIEQF MKILEEIDTL ILPENFKDSR LKRKGLVKKV QAFLAECDTV EQNICQETER LQSTNFALAE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human BAG1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.